

Debridement

The presence of a biofilm in a wound is known to be detrimental to wound healing¹⁻⁷ as is the case with the presence of necrosis or slough⁸⁻¹². Thus, the removal of these harmful aspects, debridement, is essential for creating a healing wound environment.

The type of debridement chosen depends on a number of factors, such as suitability to the patient, the type of lesion, its anatomical location, the extent of debridement required, and equipment, techniques, and clinical expertise available in a given setting.

A number of materials and techniques is available for debriding lesions and they may be combined within one patient or used sequentially. The most common and important ones are listed below.

- Surgical debridement is used in many different types of lesions and includes a wide array of materials and devices, from curettes and scalpels to highly specialized knives such as the Humby knife, the Weck dermatome and the cautery knife, which are used for large-area debridement such as excision of full thickness burns¹³⁻¹⁶. Surgical, excisional debridement is quick and safe when done by a surgical but, indeed, requires a high level of expertise¹⁷. It is also relatively expensive because it typically has to be done in an OR setting. Depending on the technique and size of the excision, blood loss may be significant as well¹⁸.
- Hydrosurgery uses high a thin, high-powered jet of saline combined with a venturi system to avoid aerosol formation, and is used in many different types of lesions¹⁹⁻²¹. It is precise and relatively quick²¹ but a certain level of expertise is required and the equipment is relatively expensive.
- Non-contact low frequency ultrasound is another way to disrupt biofilm and necrosis. It can be used in many indications although, in practice, its primary use is in skin ulcers^{22 23}. It requires a certain level of expertise and is relatively expensive but quick.
- In the USA, two types of enzymes are used for debridement. Collagenase is used by some for the debridement of superficial partial thickness burns²⁴⁻²⁶ but its use in skin ulcers is more common²⁷⁻²⁹. The compound is not very expensive and easy to use but debridement is slow with variable results³⁰ while certain topical side effects have been described^{29 31}. Bromelain, an enzyme derived from pineapples, is approved in Europe as an alternative for surgical excision in deep burns³²⁻³⁴. It has been used with favorable results. The procedures have to be done in the OR, a learning curve exists and the procedure causes pain for up to 24 hours^{32 35 36}.
- Biosurgical procedure uses a specific type of maggots. This method is very precise (since the animals only remove dead tissue) and relatively cheap and fast. Frequently, the patient has to overcome a psychological barrier, however³⁷⁻⁴⁰.
- Negative pressure wound therapy (NPWT) uses suction to remove unwanted tissue and debris from a lesion. A tight seal is created around the lesion and an external suction device is connected to the seal via a valve. Is it successfully used for both trauma and ulcers with the most extensive clinical documentation for diabetic foot ulcers⁴¹⁻⁴³ and venous leg ulcers⁴⁴⁻⁴⁶.

Although a wide range of NPWT devices exist, equipment is relatively expensive, and the seal may fail relatively frequently on certain anatomical locations.

- Surfactants, organic compounds that contain a hydrophobic and a hydrophilic region arrange themselves in spherical aggregates (micelles) when in an aqueous system. Micelles trap fragments of loose tissue and flushing the wound will remove the tissues within the micelles⁴⁷⁻⁴⁹. Surfactants are easy to use, relatively inexpensive but they take a long time for a wound to be cleaned.
- Wet-to-dry debridement uses moist gauze that, when dried out, sticks to the tissues in the wound: upon removal of the gauze the tissue is removed as well. This is one of the oldest techniques for debridement and it is easy and cheap. The major drawback is that the technique is unspecific⁵⁰⁵¹: necrosis but also viable tissue may be removed and removal of the dressing per se is typically painful⁵². These side effects make for the use of the wet-to-dry technique to be considered sub-standard practice⁵³
- Autolytic debridement relies on the bodily systems for cleaning up a wound⁵⁴. The efficacy of this method depends on the wound environment and, thus, is improved by using dressing that create a moist environment, which is shown to allow for better survivability and mobility of cells such as neutrophils and macrophages⁵⁵⁻⁵⁸. Autolytic debridement is easily promoted by using the appropriate dressings but the dressings themselves may be relatively expensive and autolysis is slow.

DEBRICHEM debrides in a way that is different from the materials and techniques mentioned above. DEBRICHEM is an active gel containing an acidic species with a potent hygroscopic action⁵⁹ which, when in contact with microorganisms and necrosis, leads to swift desiccation and carbonization of the proteins in bacteria, yeasts, and viruses, as well as of necrosis and the extracellular matrix of biofilms. The resulting dead materials rapidly slough off, thus removing elements from the wound surface that are detrimental to wound healing and allowing the development of granulation tissue and subsequent reepithelialization.

References

1. Schierle CF, De la Garza M, Mustoe TA, et al. Staphylococcal biofilms impair wound healing by delaying reepithelialization in a murine cutaneous wound model. *Wound Repair Regen* 2009;17(3):354-9. doi: 10.1111/j.1524-475X.2009.00489.x [published Online First: 2009/08/08]
2. Seth AK, Geringer MR, Gurjala AN, et al. Treatment of *Pseudomonas aeruginosa* biofilm-infected wounds with clinical wound care strategies: a quantitative study using an in vivo rabbit ear model. *Plast Reconstr Surg* 2012;129(2):262e-74e. doi: 10.1097/PRS.0b013e31823aeb3b
3. Clinton A, Carter T. Chronic Wound Biofilms: Pathogenesis and Potential Therapies. *Lab Med* 2015;46(4):277-84. doi: 10.1309/LMBNSWKUI4JPN7SO
4. Percival SL, Hill KE, Williams DW, et al. A review of the scientific evidence for biofilms in wounds. *Wound Repair Regen* 2012;20(5):647-57. doi: 10.1111/j.1524-475X.2012.00836.x [published Online First: 2012/09/19]
5. Rajpaul K. Biofilm in wound care. *Br J Community Nurs* 2015;Suppl Wound Care:S6, S8, S10-1. doi: 10.12968/bjcn.2015.20.Sup3.S6 [published Online First: 2015/03/12]
6. Schultz G, Bjarnsholt T, James GA, et al. Consensus guidelines for the identification and treatment of biofilms in chronic nonhealing wounds. *Wound Repair Regen* 2017;25(5):744-57. doi: 10.1111/wrr.12590
7. Wolcott RD, Rhoads DD, Dowd SE. Biofilms and chronic wound inflammation. *J Wound Care* 2008;17(8):333-41. doi: 10.12968/jowc.2008.17.8.30796 [published Online First: 2008/08/30]
8. Ayello EA, Cuddigan JE. Conquer chronic wounds with wound bed preparation. *Nurse Pract* 2004;29(3):8-25; quiz 26-7. doi: 10.1097/00006205-200403000-00002 [published Online First: 2004/03/17]
9. Gethin G, Cowman S, Kolbach DN. Debridement for venous leg ulcers. *Cochrane Database Syst Rev* 2015(9):CD008599. doi: 10.1002/14651858.CD008599.pub2 [published Online First: 2015/09/15]
10. Manna B, Morrison CA. Wound Debridement. *StatPearls*. Treasure Island (FL)2020.
11. Mosti G. Wound care in venous ulcers. *Phlebology* 2013;28 Suppl 1:79-85. doi: 10.1177/0268355513477015 [published Online First: 2013/03/27]
12. Turns M. Diabetic foot ulcer management: the podiatrist's perspective. *Br J Community Nurs* 2013;Suppl:S14, S16-9. doi: 10.12968/bjcn.2013.18.sup12.s14 [published Online First: 2014/05/07]
13. Evans AJ. Use of the humby knife in the excision of burns. *Br Med J* 1952;2(4791):979. doi: 10.1136/bmj.2.4791.979 [published Online First: 1952/11/01]
14. Wilder D, Rennekampff HO. [Debridement of burn wounds: rationale and options]. *Handchir Mikrochir Plast Chir* 2007;39(5):302-7. doi: 10.1055/s-2007-989227 [published Online First: 2007/11/07]
15. Hunt JL, Sato R, Baxter CR. Early tangential excision and immediate mesh autografting of deep dermal hand burns. *Ann Surg* 1979;189(2):147-51. doi: 10.1097/00000658-197902000-00004 [published Online First: 1979/02/01]
16. Mitsukawa N, Satoh K, Hosaka Y. Hemostasis by means of a cautery knife equipped with an air spray for burns over a large area. *Burns* 2006;32(6):695-7. doi: 10.1016/j.burns.2006.01.005 [published Online First: 2006/07/14]
17. Steed DL. Debridement. *Am J Surg* 2004;187(5A):71S-74S. doi: 10.1016/S0002-9610(03)00307-6 [published Online First: 2004/05/19]
18. Housinger TA, Lang D, Warden GD. A prospective study of blood loss with excisional therapy in pediatric burn patients. *J Trauma* 1993;34(2):262-3.

19. Ferrer-Sola M, Sureda-Vidal H, Altimiras-Roset J, et al. Hydrosurgery as a safe and efficient debridement method in a clinical wound unit. *J Wound Care* 2017;26(10):593-99. doi: 10.12968/jowc.2017.26.10.593 [published Online First: 2017/10/05]
20. Hong CC, Nather A, Lee JK, et al. Hydrosurgery is Effective for Debridement of Diabetic Foot Wounds. *Ann Acad Med Singap* 2014;43(8):395-9. [published Online First: 2014/09/24]
21. Madhok BM, Vowden K, Vowden P. New techniques for wound debridement. *Int Wound J* 2013;10(3):247-51. doi: 10.1111/iwj.12045 [published Online First: 2013/02/20]
22. Ennis WJ, Foremann P, Mozen N, et al. Ultrasound therapy for recalcitrant diabetic foot ulcers: results of a randomized, double-blind, controlled, multicenter study. *Ostomy Wound Manage* 2005;51(8):24-39. [published Online First: 2005/10/20]
23. Maher SF, Halverson J, Misiewicz R, et al. Low-frequency ultrasound for patients with lower leg ulcers due to chronic venous insufficiency: a report of two cases. *Ostomy Wound Manage* 2014;60(2):52-61. [published Online First: 2014/02/12]
24. Hansbrough JF, Achauer B, Dawson J, et al. Wound healing in partial-thickness burn wounds treated with collagenase ointment versus silver sulfadiazine cream. *J Burn Care Rehabil* 1995;16(3 Pt 1):241-7.
25. Ozcan C, Ergun O, Celik A, et al. Enzymatic debridement of burn wound with collagenase in children with partial-thickness burns. *Burns* 2002;28(8):791-4.
26. Patry J, Blanchette V. Enzymatic debridement with collagenase in wounds and ulcers: a systematic review and meta-analysis. *Int Wound J* 2017;14(6):1055-65. doi: 10.1111/iwj.12760 [published Online First: 2017/04/26]
27. Lee LK, Ambrus JL. Collagenase therapy for decubitus ulcers. *Geriatrics* 1975;30(5):91-3, 97-8. [published Online First: 1975/05/01]
28. Milne CT, Ciccarelli A, Lassy M. A comparison of collagenase to hydrogel dressings in maintenance debridement and wound closure. *Wounds* 2012;24(11):317-22. [published Online First: 2012/11/01]
29. Onesti MG, Fioramonti P, Fino P, et al. Effect of enzymatic debridement with two different collagenases versus mechanical debridement on chronic hard-to-heal wounds. *Int Wound J* 2016;13(6):1111-15. doi: 10.1111/iwj.12421 [published Online First: 2015/02/05]
30. Klasen HJ. A review on the nonoperative removal of necrotic tissue from burn wounds. *Burns* 2000;26(3):207-22.
31. Mosher BA, Cuddigan J, Thomas DR, et al. Outcomes of 4 methods of debridement using a decision analysis methodology. *Adv Wound Care* 1999;12(2):81-8. [published Online First: 1999/05/18]
32. Rosenberg L, Lapid O, Bogdanov-Berezovsky A, et al. Safety and efficacy of a proteolytic enzyme for enzymatic burn debridement: a preliminary report. *Burns* 2004;30(8):843-50.
33. Rosenberg L, Shoham Y, Krieger Y, et al. Minimally invasive burn care: a review of seven clinical studies of rapid and selective debridement using a bromelain-based debriding enzyme (Nexobrid(R)). *Ann Burns Fire Disasters* 2015;28(4):264-74. [published Online First: 2016/10/26]
34. Schulz A, Shoham Y, Rosenberg L, et al. Enzymatic Versus Traditional Surgical Debridement of Severely Burned Hands: A Comparison of Selectivity, Efficacy, Healing Time, and Three-Month Scar Quality. *J Burn Care Res* 2017;38(4):e745-e55. doi: 10.1097/BCR.0000000000000478 [published Online First: 2017/06/24]
35. Taussig SJ, Batkin S. Bromelain, the enzyme complex of pineapple (*Ananas comosus*) and its clinical application. An update. *J Ethnopharmacol* 1988;22(2):191-203.

36. Taussig SJ, Yokoyama MM, Chinen A, et al. Bromelain: a proteolytic enzyme and its clinical application. A review. *Hiroshima J Med Sci* 1975;24(2-3):185-93.
37. Mohd Zubir MZ, Holloway S, Mohd Noor N. Maggot Therapy in Wound Healing: A Systematic Review. *Int J Environ Res Public Health* 2020;17(17) doi: 10.3390/ijerph17176103 [published Online First: 2020/08/23]
38. Munro S, Hadid A, Rahmani MJH. Maggots in the management of ulcer care. *BMJ Case Rep* 2017;2017 doi: 10.1136/bcr-2017-220462 [published Online First: 2017/05/28]
39. Webb R. For centuries in wound healing. *J Wound Care* 2017;26(3):77. doi: 10.12968/jowc.2017.26.3.77 [published Online First: 2017/03/10]
40. Falch BM, de Weerd L, Sundsfjord A. [Maggot therapy in wound management]. *Tidsskr Nor Laegeforen* 2009;129(18):1864-7. doi: 10.4045/tidsskr.08.0656 [published Online First: 2009/10/22]
41. Alavi A, Sibbald RG, Mayer D, et al. Diabetic foot ulcers: Part II. Management. *J Am Acad Dermatol* 2014;70(1):21 e1-24; quiz 45-6. doi: 10.1016/j.jaad.2013.07.048 [published Online First: 2013/12/21]
42. Andros G, Armstrong DG, Attinger CE, et al. Consensus statement on negative pressure wound therapy (V.A.C. Therapy) for the management of diabetic foot wounds. *Ostomy Wound Manage* 2006;Suppl:1-32. [published Online First: 2006/09/30]
43. Blume PA, Walters J, Payne W, et al. Comparison of negative pressure wound therapy using vacuum-assisted closure with advanced moist wound therapy in the treatment of diabetic foot ulcers: a multicenter randomized controlled trial. *Diabetes Care* 2008;31(4):631-6. doi: 10.2337/dc07-2196 [published Online First: 2007/12/29]
44. Karavan M, Olerud J, Bouldin E, et al. Evidence-based chronic ulcer care and lower limb outcomes among Pacific Northwest veterans. *Wound Repair Regen* 2015;23(5):745-52. doi: 10.1111/wrr.12341 [published Online First: 2015/07/15]
45. Richmond NA, Maderal AD, Vivas AC. Evidence-based management of common chronic lower extremity ulcers. *Dermatol Ther* 2013;26(3):187-96. doi: 10.1111/dth.12051 [published Online First: 2013/06/08]
46. Ross RE, Aflaki P, Gendics C, et al. Complex lower extremity wounds treated with skin grafts and NPWT: a retrospective review. *J Wound Care* 2011;20(10):490, 92-5. doi: 10.12968/jowc.2011.20.10.490 [published Online First: 2011/11/10]
47. Pittinger TP, Curran D, Hermans MH. The treatment of paediatric burns with concentrated surfactant gel technology: a case series. *J Wound Care* 2020;29(Sup6):S12-S17. doi: 10.12968/jowc.2020.29.Sup6.S12 [published Online First: 2020/06/13]
48. Powers JG, Higham C, Broussard K, et al. Wound healing and treating wounds: Chronic wound care and management. *J Am Acad Dermatol* 2016;74(4):607-25; quiz 25-6. doi: 10.1016/j.jaad.2015.08.070 [published Online First: 2016/03/17]
49. Woo K, Hill R, LeBlanc K, et al. Effect of a surfactant-based gel on patient quality of life. *J Wound Care* 2018;27(10):664-78. doi: 10.12968/jowc.2018.27.10.664 [published Online First: 2018/10/18]
50. Halim AS, Khoo TL, Saad AZ. Wound bed preparation from a clinical perspective. *Indian J Plast Surg* 2012;45(2):193-202. doi: 10.4103/0970-0358.101277 [published Online First: 2012/11/2]
51. Kammerlander G, Andriessen A, Asmussen P, et al. Role of the wet-to-dry phase of cleansing in preparing the chronic wound bed for dressing application. *J Wound Care* 2005;14(8):349-52. doi: 10.12968/jowc.2005.14.8.26824 [published Online First: 2005/09/24]

52. Cowan LJ, Stechmiller J. Prevalence of wet-to-dry dressings in wound care. *Adv Skin Wound Care* 2009;22(12):567-73. doi: 10.1097/01.ASW.0000363469.25740.74 [published Online First: 2009/11/26]
53. Wodash AJ. Wet-to-Dry Dressings Do Not Provide Moist Wound Healing. *J Am Coll Clin Wound Spec* 2012;4(3):63-6. doi: 10.1016/j.jccw.2013.08.001 [published Online First: 2012/09/01]
54. Mancini S, Cuomo R, Poggialini M, et al. Autolytic debridement and management of bacterial load with an occlusive hydroactive dressing impregnated with polyhexamethylene biguanide. *Acta Biomed* 2018;88(4):409-13. doi: 10.23750/abm.v88i4.5802 [published Online First: 2018/01/20]
55. Hermans MH. HydroColloid dressing (Duoderm) for the treatment of superficial and deep partial thickness burns. *Scand J Plast Reconstr Surg Hand Surg* 1987;21(3):283-5. doi: 10.3109/02844318709086461 [published Online First: 1987/01/01]
56. Hermans MH. Hydrocolloid dressing versus tulle gauze in the treatment of abrasions in cyclists. *Int J Sports Med* 1991;12(6):581-4. doi: 10.1055/s-2007-1024738 [published Online First: 1991/12/01]
57. McFarland A, Smith F. Wound debridement: a clinical update. *Nurs Stand* 2014;28(52):51-8. doi: 10.7748/ns.28.52.51.e9171 [published Online First: 2014/08/28]
58. Reyzelman AM, Vartivarian M. Evidence of Intensive Autolytic Debridement With a Self-Adaptive Wound Dressing. *Wounds* 2015;27(8):229-35. [published Online First: 2015/08/19]
59. Compositions for removing necrotic or infected tissues from body lesions and from oral cavity. Patent Cooperation Treaty;Nr. PCT/IB2019/051146:Patent Pending.